

# UNLEASH YOUR ULTIMATE SEQUENCING SPEED

## Benchtop Genetic Sequencer **DNBSEQ-G99**



- **Rapid sequencing**  
Only 12 hrs for PE150 (from loading to FASTQ).
- **Flexible throughput**  
Independent loading and running of dual flow cells.
- **Bioinformatics integrated**  
Option to include built-in bioinformatics module to support sequencing and advanced analysis in a single machine.



MGI Tech Co., Ltd.

Building 11, Beishan Industrial Zone,  
Yantian District, Shenzhen, CHINA  
518083

Version: August 2023 | MGPA0603002

Information in this brochure is updated to [08/02/2023] and only for your reference. In no event shall the brochure be regarded as warranty or commitment made by MGI Tech Co., Ltd. All rights and obligations shall be subject to the final executed agreement.




[en.mgi-tech.com](http://en.mgi-tech.com)

[MGI-service@mgi-tech.com](mailto:MGI-service@mgi-tech.com)

\*Unless otherwise informed, StandardMPS and CoolMPS sequencing reagents, and sequencers for use with such reagents are not available in Germany, Spain, UK, Sweden, Italy, Czech Republic, Switzerland and Hong Kong (CoolMPS is available in Hong Kong).

MGI' S PROPRIETARY  
**DNBSEQ™**  
TECHNOLOGY

# Genetic Sequencer **DNBSEQ-G99**



DNBSEQ-G99 is developed based on MGI' s core DNBSEQ™ sequencing technology. Enabled by innovations in biochemistry, optics, fluidics, temperature control, and other core systems, DNBSEQ-G99 boasts the fastest speed amongst all medium-to-low throughput sequencers globally. DNBSEQ-G99 is especially applicable for targeted oncology panel sequencing, infectious disease sequencing, oncology methylation sequencing, small whole-genome sequencing, low-depth whole genome sequencing, individual identification, 16s metagenomics sequencing, small panel sequencing of 24-28 samples, or whole-exome sequencing of 1-4 samples.

Powered by 4-color sequencing technology, DNBSEQ-G99 also comes with an optional built-in bioinformatics module, which allows advanced analysis to begin automatically after the sequencing run. This facilitates a tremendously efficient and simple workflow, thus accelerating the application of omics technology to advance global life sciences and clinical research.


# Designed for Simplicity

Newly designed flow cell, reagent cartridge, and user interface are introduced in DNBSEQ-G99, providing laboratory personnel with unparalleled ease and peace of mind in the entire sequencing workflow. A built-in bioinformatics module can also be included, achieving from sample to report all in one equipment.



## Novel Reagent Cartridge Design

- Pre-loaded reagents within cartridge, one-step operation: simply press to load
- Sequencing and cleaning cartridges combined 2-in-1, cleaning initiates automatically after run.



## Intelligent Interaction

- Visualize the entire sequencing process in real-time
- Intuitive animations are included to guide flow cell loading, minimizing operational errors



## Built-in Bioinformatics Module

- Advanced analysis begins automatically after run, and supports Bioanalysis by Sequencing (BBS) mode
- ZLIMS-compatible, achieve efficient workflow management and local data output



## Data Security

- Designed based on GDPR privacy protection requirement
- Secure storage to safeguard your sensitive data




# Flexible

## Customize your run

DNBSEQ-G99 is the only medium-to-low throughput sequencer with a dual flow cell loading configuration, providing ultimate flexibility in sequencing throughput. The dual flow cells can be operated with different read lengths independently or concurrently. Three flow cell loading modes are supported on DNBSEQ-G99: Single, Dual concurrent, or Dual independent. Laboratory technicians can decide number of flow cells to operate in accordance with the sample size and requirements.



**Dual flow cells**  
Plug-and-play as needed



**Reagent cartridge**  
Two independent cartridges for two flow cells




**Bioinformatics module**  
Analysis automatically begins after each flow cell run

# Speed


for your time-sensitive needs

To achieve the fastest sequencing speed in its range, DNBSEQ-G99 incorporates innovatively optimized flow cell, biochemistry process, fluidics, optics, temperature control and other core systems. The result is supreme sequencing efficiency, while data quality remains top-notch.




## High Density Flow Cell

- 600nm pitch high density patterning
- 68% more DNB loading per unit area
- Novel triangular configuration




## Super Fast Biochemistry

- 10s rapid fluorescence reaction
- Biochemistry incubation reaction sped up from minutes-range to seconds-range



## Surpassing the Optical Diffraction Limit

- In-house developed ultra-high quality objective lens
- Improvement of signal capture efficiency by decreasing scan area



## Rapid Temperature Control

- ~7 °C/s for heating and cooling
- Doubles the heating and cooling speed of conventional PCR instrument



# Superior Performance


DNBSEQ-G99 delivers uncompromised high data quality. A multitude of applications can be executed on DNBSEQ-G99, such as targeted sequencing, small genome and Low pass WGS sequencing, etc.

In addition, DNBSEQ-G99A supports the retrieval of data at intermediate time points under the Bioanalysis by Sequencing (BBS) mode. Users can obtain the first batch of summary report as quick as 2.5 hrs from the start of sequencing run (read length: SE40).

Method	Application	Recommended read length	Data size per Sample	Samples per Run
Targeted Capture/ Multiplex PCR	Oncology panel	PE100, PE150	Small panel: ~1 Gb/sample	24/FC, 48/RUN
	Hereditary disease small panel (Thalassemia, deafness, etc.)	PE150	Deafness: ~5 Gb/sample Thalassemia: ~0.2 M reads/sample	Deafness: 4/FC, 8/RUN Thalassemia: 400/FC, 800/RUN
	ATOPIlex panel (respiratory disease, SARS-CoV-2, etc.)	PE100, PE150	Respiratory tract panel: 5 M reads/sample COVID-19 panel: 5 M reads/sample	16/FC, 32/RUN
	WES	PE150	~15 Gb/sample	1-2/FC, 2-4/RUN
Methylation Analysis	Oncology targeted methylation panel	PE150	~5 Gb/sample	4/FC, 8/RUN
Small Genome Sequencing	Metagenomics for pathogen detection	SE50, SE100	Meta: 20 M reads/sample	4/FC, 8/RUN
	Microbial WGS	PE100, PE150	Isolated bacteria: ~1 Gb/sample	16-24/FC, 32-48/RUN
	16s V3-V4 sequencing	PE300	≥0.1 M reads/sample	576/FC, 1152/RUN
Low pass whole-genome sequencing	NIPT	SE50	NIPT/PGS: ~10 M reads/sample	8/FC, 16/RUN
	PGS	SE50		
RNA sequencing	Expression profiling Transcriptome	SE50	Expression profiling: ~25 M reads/sample	Expression profiling: 3/FC, 6/RUN
		PE150	Transcriptome: ~6 Gb/sample	Transcriptome: 4/FC, 8/RUN

① Recommended data output and sample numbers are only for reference, actual application will require optimisation adjustments.

② SE400 will be released in 2023.

③  Recommend method

## Performance Parameters

Maximum number of Flow cells	Lanes/ Flow cell	Effective Reads*/ Flow Cell	Supported Reads Lengths <sup>+</sup>	Data Output	Q30 <sup>**</sup>	Run Time
2	1	80M	SE100/PE50	8~16G	>90%	5h
			PE150	24~48G	>85%	12h
			APP-C SE100 <sup>***</sup>	8~16G	>90%	5h
			APP-C PE150	24~48G	>85%	12h
			PE300	48~96G	>85%	30h

\* The effective reads are based on the sequencing of an internal standard library. Actual output may vary depending on sample type and library preparation method.

<sup>+</sup> The percentage of bases above Q30 is the average of an internal standard library over the entire run. Actual performance is affected by factors such as sample type, library quality, and insert fragment length.

<sup>\*\*\*</sup> Estimated release in 2023 Q1

DNBSEQ-G99 also supports SE50 and PE100 sequencing, and the existing kits can support SE50, PE100 read length sequencing.

## Available Models



# Oncology Application

## Low Frequency Variants Detection

### Experiment Scheme

**Sample:** Lung cancer ctDNA standards, diluted to 1%, 0.5%, 0.2%, and 0.1% variant ctDNA samples

**Library prep:** Targeted capture kit from third party

**Sequencing strategy:** PE100 dual-barcode sequencing, 4 repeated runs were tested

**Objective:** To test DNBSEQ-G99 variant detection capability

### Sequencing Summary

The 4 runs generated 126 M reads on average, Q30 >93%, with excellent uniformity observed.

### Analysis Summary





- 100% detection of SNV mutation sites in the samples (1%, 0.5%, 0.2%, 0.1%).



### Sequencing Result

	Total reads (M)	Q30	EstErr(%)	Time (h)
Mean value	126.22	94.00	0.22	9.15
Standard deviation	7.41	0.52	0.02	0.11

### Analysis Result



# Oncology Application

## Methylation Sequencing (Targeted Capture)

### Experiment Scheme

**Sample:** Fragmented DNA from 4 cell lines

**Library prep:** Targeted capture double stranded library kit from third party

**Sequencing strategy:** PE100 dual-barcode sequencing, 2 repeated runs were tested

**Objective:** To evaluate DNBSEQ-G99's compatability with low-diversity libraries, and data uniformity for targeted methylation regions.

### Sequencing Summary

Without addition of a spike-in balanced library, 2 runs had an average output of 101 M reads, with Q30 >86 %, showing good compatibility with the low diversity libraries.



### Analysis Summary

- Highly concordant average methylation fractions (AMF) measured between both runs for all tested samples.
- For results of the same samples on different instruments: AMF of the samples is highly consistent.

### Sequencing Results

	Total reads (M)	Q30 (%)	SplitRate (%)	Time (h)
Run 1	103.12	86.42	98.36	9.20
Run 2	99.08	90.57	98.02	9.15

### Analysis Results





# Small Genome Sequencing

## Pathogen detection

### Experiment Scheme

**Sample:** 4-pooled 1% reference microbial community standards  
**Library prep:** MGIEasy FS DNA Library Prep Set  
**Sequencing strategy:** PE100 single-barcode sequencing  
**Objective:** Assess DNBSEQ-G99's capability to identify unknown pathogens

### Sequencing Summary

Output of 110M reads, Q30 >95 %, exceeding data amount required for analysis.

### Analysis Summary

- Pathogen fast identification (PFI) was used for analysis, and the pathogen identification results were consistent with reference microbial community in terms of detected species and abundance.
- Fluctuation in abundance CV was lower than 2 %, indicating high accuracy.

### Sequencing Results

	Total reads (M)	Q30 (%)	SplitRate (%)	Time (h)
Output	110.17	95.22	96.95	9

### Analysis Results

Species	Sample 1	Sample 2	Sample 3	Sample 4	Standard abundance	mean	SD	CV
<i>Salmonella enterica</i>	15.87%	15.60%	15.72%	15.87%	12.00%	15.77%	0.13%	0.82%
<i>Pseudomonas aeruginosa</i>	14.12%	14.00%	13.83%	14.12%	12.00%	14.02%	0.14%	1.00%
<i>Bacillus subtilis</i>	13.30%	13.14%	13.46%	13.30%	12.00%	13.30%	0.13%	0.98%
<i>Escherichia coli</i>	11.75%	12.32%	11.90%	11.75%	12.00%	11.93%	0.27%	2.26%
<i>Enterococcus faecalis</i>	11.13%	11.30%	11.24%	11.13%	12.00%	11.20%	0.08%	0.71%
<i>Listeria monocytogenes</i>	11.12%	11.02%	11.12%	11.12%	12.00%	11.10%	0.05%	0.45%
<i>Staphylococcus aureus</i>	10.18%	10.23%	10.30%	10.18%	12.00%	10.22%	0.06%	0.59%
<i>Limosilactobacillus fermentum</i>	9.52%	9.48%	9.46%	9.52%	12.00%	9.50%	0.03%	0.32%
<i>Cryptococcus neoformans</i>	1.49%	1.49%	1.52%	1.49%	2.00%	1.50%	0.02%	1.33%
<i>Saccharomyces cerevisiae</i>	1.47%	1.42%	1.45%	1.47%	2.00%	1.45%	0.02%	1.38%

# Small Genome Sequencing

## Phage Assembly

### Experiment Scheme

**Sample:** 16 pure bacteriophage cultures

**Library prep:** MGIEasy universal DNA library prep set

**Test Strategy:** PE150 dual-barcode

**Test Purpose:** Assess DNBSEQ-G99 capability for assembling whole bacteriophage genome

### Sequencing Summary

Output of 106M reads, Q30 >92 %, exceeding data amount required for analysis.


### Analysis Summary

- Conserved protein genes for phage capsid, portal, and terminase were successfully detected in the DNBSEQ-G99 assemblies.
- DNBSEQ-G99 outperforms other platform for bacteriophage whole genome assembly, achieving higher assembly integrity.

### Sequencing Results

	Total reads (M)	Q30 (%)	SplitRate (%)	Time (h)
Output	106.99	92.32	97.97	12

### Analysis Results



# Small Genome Sequencing

## 16s Sequencing

### Experiment Scheme

**Sample:** Zymobiomics D6305 reference standard samples  
**Library prep:** ATOplex 16S V3 V4 rDNA Library Preparation Set  
**Sequencing strategy:** PE300 dual barcode sequencing  
**Objective:** To evaluate the data quality of DNBSEQ-G99 for 16s samples

### Sequencing Summary

4 runs had an average output of 93.45 M, with Q30>90%, which can meet the demand of bioinformatics analysis.


### Analysis Summary

- The OUT results showed that the abundance at the genus level was highly consistent with the expected results.
- For the same sample in different runs, the correlation coefficient  $R^2 > 0.99$  between parallel sequencing libraries, which showed high consistency of test results.

### Sequencing Result

	Total reads (M)	Q30 (%)	SplitRate (%)
Run-1	89.89	89.93	96.28
Run-2	94.97	90.21	96.65
Run-3	99.68	90.2	96.8
Run-4	89.25	89.82	95.34
Average	93.45	90.04	96.27

### Analysis Result





## Hardware Specifications

<b>Model</b>	DNBSEQ-G99 DNBSEQ-G99A	Outputs FASTQ files Equipped with bioinformatics module for advanced analysis
<b>Dimensions/Net Weight</b>	607*680*640 mm/~140 kg	
<b>Power</b>	Rated Voltage Rated frequency Rated Power	100 V-240 V 50/60 Hz 1000 VA, [working current]: ≥10 A
<b>Touch Screen</b>	LCD touch screen	
	Touch screen size Touch screen resolution	21.5 inch 1920 x 1080
<b>Maximum Sound Pressure</b>	75 dB(A)	
<b>Shell Protection Grade</b>	IPX0	
<b>Operating Environment Requirements</b>	Temperature Relative Humidity Atmospheric Pressure Maximum Altitude (above sea level)	19-30 °C 20-80 %RH 70 kPa-106 kPa 3000 m
<b>Computer Configurations</b>	CPU Internal Storage HDD Operating System	Intel i9-10900e 2.80 GHz 64 GB 6 TB Windows 10
<b>Bioinformatics Module Configurations</b>	CPU Memory System Disk Cache Disk Storage Disk Ethernet	Intel Xeon 5220S 18C/36T 2.7GHz * 2 256 GB 960 GB 960 GB 32 TB Gigabit Ethernet RJ45 * 2

\* The maximum sound pressure is measured and calculated at any position with the maximum sound pressure level 1m away from the housing during normal use

\*\* For indoor use only

\*\*\* Support computer configuration and system version upgrade

## Ordering Information

RUO\*

Cat. No	Product Name
900-000607-00	DNBSEQ-G99RS
900-000609-00	DNBSEQ-G99ARS
940-000409-00	High-throughput Sequencing Set (G99 SM FCL SE100/PE50)
940-000410-00	High-throughput Sequencing Set (G99 SM FCL PE150)
940-000415-00	High-throughput Sequencing Set (G99 SM FCL PE300)
940-000413-00	High-throughput Sequencing Set (G99 SM FCL APP-C PE150)
940-000520-00	High-throughput Sequencing Set (G99 SM App-C FCL SE100)
940-000624-00	DNBSEQ-G99 Cleaning Reagent Kit
Selected as needed	UPS

\* For research use only. Not for use in diagnostic procedures

IVD

Cat. No	Product Name
900-000612-00	DNBSEQ-G99
900-000628-00	DNBSEQ-G99A
940-000428-00	Universal Sequencing Reaction Kit (G99 SM FCL SE100/PE50)
940-000431-00	Universal Sequencing Reaction Kit (G99 SM FCL PE150)
940-000434-00	Universal Sequencing Reaction Kit (G99 SM App-C FCL PE150)
940-000525-00	Universal Sequencing Reaction Kit (G99 SM App-C FCL SE100)

## MGI Global Presence

### Technical Support Available Globally

The MGI technical support team has a complete global coverage, including technical service centers in major international regions and multiple locations to maximize customer satisfaction.



Multiple local technical support centers and Customer Experience Centers (CEC) worldwide are established to provide timely and effective technical support and training.



Multiple local spare part centers worldwide are established to ensure a sufficient supply of parts for machine maintenance.



Online technical support accessible worldwide with a fully functioning call center (Toll-Free Hotline 4000-688-114) (9:00 AM-12:00 PM, 13:00 PM-18:00 PM, Beijing time (GMT+8), workday).

### Comprehensive Instrument Service and Warranty Plans Globally



Multiple local warehouses worldwide are established to ensure a timely and effective supply of maintenance parts.



Providing installation and system verification services (inclusive of necessary reagents and consumables).



Responsible for any manufacturing defects or faults on the system within the warranty.



Providing instrument preventive maintenance within warranty, along with a variety of available extended warranty support plans.

## MGI Genetic Sequencers



### DNBSEQ-E25

Reads per flow cell: 25M  
Number of flow cells: 1  
Data output: 2.5-7.5Gb



### DNBSEQ-G99

Reads per flow cell: 80M  
Number of flow cells: 2  
Data output: 8-96GB



### DNBSEQ-G50

Reads per flow cell: 100-500M  
Number of flow cells: 1  
Data output: 10-150GB



### DNBSEQ-G400

Reads per flow cell: 300-1800M  
Number of flow cells: 2  
Data output: 55-1440GB



### DNBSEQ-T7

Reads per flow cell: 5800M  
Number of flow cells: 4  
Data output: 1-7Tb



### DNBSEQ-T20x2

Reads: 35-40Billion  
Number of sides: 6  
Data output: 42-72TB



# About

## MGI Tech Co., Ltd.

**2,800+**  
Employees

MGI Tech Co., Ltd. (referred to as MGI) is committed to building core tools and technology to lead life science through intelligent innovation. With a focus on R&D, production and sales of DNA sequencing instruments, reagents, and related products, MGI provides real-time, panoramic, and life course equipment and systems for precision medicine, precision agriculture, precision healthcare and other relevant industries. MGI is a leading producer of clinical high-throughput gene sequencers, and its multi-omics platforms include genetic sequencing, medical imaging, and laboratory automation.

**36%**  
R&D Personnel

As of December 31, 2022, MGI has more than 2,800 employees, and 36% of whom are R&D personnel. Founded in 2016, MGI operates in more than 90 countries and regions, serving more than 2,000 customers. It has established scientific research and production bases, global training and service network in many countries and regions around the world. MGI is one of two companies in the world that can independently develop and mass-produce low-, medium- and high-throughput clinical gene sequencers from GB to TB. Providing real-time, comprehensive, life course solutions, its vision is to lead life science innovation.

**2,000+**  
Customers

**Vision**  
Leading Life Science Innovation

**90+**  
Countries & Regions

**Mission**  
To Develop and Promote Advanced Life Science Tools for Future Healthcare