GLDH VET test kit

For veterinary use only!

Veterinary test kit for quantitative in vitro determination of GLDH in serum or lithium heparin plasma with a solo analyser.

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Test kit preparation: allow single test at least 10 minutes to warm up to room temperature (20-25 °C) by placing the test into the test kit rack. Return the test kit package to the refrigerator.

Summary

Glutamatedehydrogenase (GLDH) is a liver specific enzyme which is located in mitochondria of centrolobular hepatocytes. Due to this localization, GLDH dominantly increases in hypoxic states of the liver as well as in pronounced necrosis of hepatocytes.

Increased release of GLDH from hepatocytes can further be observed in i.e. inflammatory, neoplastic, or infectious liver pathology as well as in all diseases leading to liver cell degeneration.

Method

Glutamatedehydrogenase catalyzes the reduction of Oxolutarat and at the same time the oxidation of NADH to NAD.

The reaction speed of the absorption increase is directly proportional to the GLDH activity which is measured on the SOLO system at 340 nm.

Measurement Range

6 - 120 U/L (Lot dependent)

Sample Material

Use 100µl lithium heparin plasma or serum samples. (Avoid lipaemic samples)

In case samples need to be kept, store at a dark place only!

Sample material stability:

- Stored at 20 °C: 4 hours
- Stored at 2-8 °C: 1 day
- Stored at -20 °C: 4 weeks

Test Kit

ERS Cuvette filled with: Trietanolamin, Alpha Ketoglutarat, Ammoniumacetat EDTA, ADP, LDH.

ERS Cap filled with: NADH

Stability and Storage

Stable until the expiration date stated on the label when stored in unopened vacuum package at 2-8 °C. Opening the vacuum package may limit the reagent stability to three months (stored at 2-8 °C) from the date of opening.
DO NOT FREEZE!

Warnings and Precautions

DO NOT INGEST! Avoid contact with skin and eyes. Observe all necessary precautions for the use of laboratory reagents.

Waste Management

Always refer to local legal requirements.

Reference Range

Dog:	< 13 U/L
Cat:	< 10 U/L
Rabbit:	< 19 U/L
Guinea pig:	< 17 U/L
Horse:	< 13 U/L
Bovine:	< 15 U/L
Sheep:	< 11 U/L
Pig:	< 7 U/L
Goat:	< 20 U/L

It is recommended that each laboratory establishes its own reference ranges.

Quality Control

For internal quality control the Eurolyser GLDH VET control kit is recommended. Order number: VT 2000

Precision

Reproducibility "day to day":
Control; N=20; mean = 24.4 U/L; CV = 8.02%;

Correlation

y (Solo GLDH) = 1.0414 x (COBAS GLDH) + 4.2329;
R² = 0.9494.

Interferences

The test system has been analysed for various interferences. Criterion was the recovery within 10% of initial values.

Ascorbic Acid	30 mg/dl
Bilirubin	60 mg/dl
Hemoglobin	500 mg/dl
Triglyceride	400 mg/dl

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English

ATTENTION!

Allow single test at least 10 minutes to warm up to room temperature (20-25 °C)!

1. Preparing the test system

- 1.1 Insert RFID card
- 1.2 Place ERS cuvette into test kit rack
- 1.3 Place ERS cap into test kit rack
- 1.4 Press measurement button and enter all required information using the touch screen

2. Sample preparation

- 2.1 Pipette 100µl sample material from centrifuged sample tube

3. Test processing

- 3.1 Dispense 100µl sample INTO THE LIQUID in the ERS cuvette
- 3.2 Firmly apply ERS cap on ERS cuvette
- 3.3 Place ERS cartridge into the solo analyser
- 3.4 Start automatic sample processing by pressing the  start button on the solo analyser