cCRP VET test kit

For veterinary use only!

Veterinary test kit for quantitative in vitro determination of canine CRP (cCRP) in serum, EDTA plasma or lithium heparin plasma with a solo analyser.

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Kit size

16 tests/kit
6 tests/kit
1 x 0.5 ml (high level)

Test kit preparation: allow single test at least 10 minutes to warm up to room temperature (> 20 °C) by placing the test into the test kit rack. Return the test kit package to the refrigerator.

Summary

Veterinary in vitro test kit designed to measure canine C-reactive protein (cCRP) in serum and plasma. C-reactive protein (CRP) is an acute phase protein. The serum level of CRP rises in animals that suffer from different types of infection or inflammation. It correlates with degree and activity of the causative disease. Measurements of serum CRP are used to aid the evaluation of the inflammation associated with infection, tissues damage and other clinical disorders. Increased levels of CRP were described in dogs that suffered from pneumonia, pyometra, pancreatitis, parvovirus infection or traumatism. Dogs with an increased number of leukocytes had significantly higher CRP levels than dogs with noninflammatory leukograms. Elevated CRP levels were also found in dogs with arthritis, thrombophlebitis or proctitis and in dogs infected with Bordetella bronchiseptica or Ehrlichia canis.

Method

Immunoturbidimetric determination at 546nm of the sample turbidity which is proportional to the concentration of cCRP present in the original specimen.

Measurement Range

10 - 200 mg/l (Lot dependent)

Samples with concentrations higher than the upper limit of the measurement range must be diluted 1:4 with physiological saline (0.9 % NaCl solution), e.g. 20 µl sample + 80 µl 0.9 % NaCl solution, and the result multiplied by 5.

Values up to 800 mg/l show no prozone effect.

Sample Material

Use 5µl serum, lithium heparin plasma or EDTA plasma. Sample stability testing showed that cCRP (in serum) was stable for two weeks when stored at 2-8 °C.

Test Kit

ERS Cuvette filled with reaction buffer.
ERS Cap filled with immunoparticles.

Stability and Storage

Stable until the expiration date stated on the label when stored in unopened vacuum package at 2-8 °C.
Opening the vacuum package may limit the reagent stability to three months (stored at 2-8 °C) from the date of opening.
DO NOT FREEZE!

Warnings and Precautions

DO NOT INGEST! Avoid contact with skin and eyes. Observe all necessary precautions for the use of laboratory reagents.

Waste Management

Always refer to local legal requirements.

Reference Range

Dog: < 10 mg/l (< 1 mg/dl)

It is recommended that each laboratory establishes its own reference range.

Quality Control

For internal quality control the cCRP VET control kit is recommended. Order number: VT 1000

Precision

Reproducibility within run:
Dog serum; N = 20; mean = 77 mg/l; CV = 4.7 %;

Correlation

Correlation with Gentian cCRP, 42 samples splitted in serum and lithium heparin:
 y (Eurolyser) = 1.019 x (Gentian cCRP) - 0.76; R^2 = 0.98.

Interferences

The test system has been analysed for various interferences. Criterion was the recovery within 10% of initial values.

Hemoglobin 500 mg/dl
Intralipid 1000 mg/dl

References

- Hillström, Hagman, Söder, Haggström, Ljungvall & Kjelgaard-Hansen 2015, Validation and application of a canine-specific automated high-sensitivity C-reactive protein assay

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Bezeichnung

cCRP VET test kit
cCRP VET test kit
cCRP VET control kit

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16 Tests/Packung
6 Tests/Packung
1 x 0,5 ml (high level)

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
Hämoglobin 500 mg/dl
Intralipid 1000 mg/dl

Referenzen

- Hillström, Hagman, Söder, Haggström, Ljungvall & Kjølgaard-Hansen 2015, Validation and application of a canine-specific automated high-sensitivity C-reactive protein assay

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





English

ATTENTION!

Allow single test at least 10 minutes to warm up to room temperature (>20 °C)!

1. Preparing the test system

- 1.1 Insert RFID card
- 1.2 Place ERS cuvette into test kit rack
- 1.3 Place ERS cap into test kit rack
- 1.4 Press measurement button and enter required information using the touch screen

2. Sample preparation

- 2.1 Pipette 5µl sample material from centrifuged sample tube

3. Test processing

- 3.1 Dispense 5µl sample INTO THE LIQUID in the ERS cuvette
- 3.2 Firmly apply ERS cap on ERS cuvette
- 3.3 Place ERS cartridge into the solo analyser
- 3.4 Start automatic sample processing by pressing the  start button on the solo analyser