# Fructosamin (gold standard)

## VET test kit

### For veterinary use only!

Veterinary test kit for quantitative in vitro determination of fructosamin in serum, EDTA plasma or lithium heparin plasma with a solo analyser.

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### Kit size

16 tests/kit  
6 tests/kit  
2 x 1 ml (low/high)

**Test kit preparation: allow single test at least 10 minutes to warm up to room temperature (> 20 °C) by placing the test into the test kit rack. Return the test kit package to the refrigerator.**

## Summary

Patients with diabetes mellitus may be defined as having a disturbance in the oxidation and use of glucose for energy production. This may be secondary to decreased production or release of insulin by beta cells in the pancreatic islets of Langerhans, or due to an abnormal tissue response to the presence of insulin. This disease is relatively common in both dogs and cats and is treated with insulin injections and dietary restriction. Regulation of insulin administration is important in controlling diabetes. Fructosamine is formed by the glycosylation of circulating proteins and is related to circulating glucose concentration. Fructosamine measurement can aid in the diagnosis of diabetes mellitus and in the assessment of blood glucose control in treated diabetic patients.

Eurolyser's "Gold standard Fructosamine" assay is an enzymatic assay which is much more specific for glycated serum protein (Fructosamine) than the old/traditional NBT based Fructosamine assay which may be interfered by endogenous reducing substances. The traditional method (NBT) was found to be less specific for fructosamine. NBT reacts with various endogenous substances of reducing activity such as thiol groups, NADH and ascorbate. Studies showed that only about half of the reducing activity (Fructosamine) was due to specific non enzymatic glycation of proteins, and the remaining unspecific activity varied from serum to serum. Therefore the novel Eurolyser "Gold standard Fructosamine" is the assay of the first choice offering improved specificity and reliability as it does not suffer from nonspecific interferences.

## Method

Photometric measurement of fructosamine at 546nm as enzymatic kinetic test.

## Measurement Range

0 - 2600 µmol/l (Lot dependent)

## Sample Material

Use 20µl fresh EDTA plasma, lithium heparin plasma or serum samples.

Samples can be stored at 2-8 °C for two weeks or at -20 °C for up to four weeks.

## Test Kit

ERS Cuvette filled with: Enzyme/substrate reagent containing 4-AA, Proteinase K and stabilizers

ERS Cap filled with: Enzyme/substrate reagent containing Fructosaminase, TOOS, HRP, Geneticin

## Stability and Storage

Stable until the expiration date stated on the label when stored in unopened vacuum package at 2-8 °C.

Opening the vacuum package may limit the reagent stability to three months (stored at 2-8 °C) from the date of opening.

DO NOT FREEZE!

## Warnings and Precautions

DO NOT INGEST! Avoid contact with skin and eyes. Observe all necessary precautions for the use of laboratory reagents.

## Waste Management

Always refer to local legal requirements.

## Reference Range

Cat:	160 - 365 µmol/l
Dog:	180 - 365 µmol/l
Mouse:	160 - 365 µmol/l
Rodent:	160 - 365 µmol/l
Horse:	175 - 280 µmol/l

It is recommended that each laboratory establishes its own reference ranges.

## Quality Control

For internal quality control the Eurolyser Fructosamin VET control kit is recommended. Order number: VT 1100

## Precision

Reproducibility within run:

Control; N = 80; mean = 213.26 µmol/l; CV = 3.6 %;

## Correlation

Dog sample correlation:

y (Solo Fruct) = 0.9633 x(PENTRA400) + 0.0899; R<sup>2</sup> = 0.976.

Cat sample correlation:

y (Solo Fruct) = 1.0414 x(PENTRA400) - 15.061; R<sup>2</sup> = 0.9753.

## Interferences

The test system has been analysed for various interferences. Criterion was the recovery within 10% of initial values.

Ascorbic Acid	5 mg/dl
Bilirubin	7.5 mg/dl
Bilirubin Conjugated	5 mg/dl
Glucose	2400 mg/dl
Hemoglobin	200 mg/dl
Triglyceride	2000 mg/dl
Uric Acid	35 mg/dl

## References

1. Armbruster DA, Fructosamine: Structure, Analysis and Clinical Usefulness. Clin. Chem. 1987; 33 (12): 2153-2163.
2. Kouzuma, T. et al. An Enzymatic Method for the Measurement of Glycated Albumin in Biological Samples. Clin. Chim. Acta 2002; 324: 61-71

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
### Referenzen


1. Ambruster DA, Fructosamine: Structure, Analysis and Clinical Usefulness. Clin. Chem. 1987; 33 (12): 2153-2163.
2. Kouzuma, T. et al. An Enzymatic Method for the Measurement of Glycated Albumin in Biological Samples. Clin. Chim. Acta 2002; 324: 61-71

# Durchführung eines Fructosamin (gold standard) VET Tests Performing a Fructosamin (gold standard) VET test

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






## English

### ATTENTION!

Allow single test at least 10 minutes to warm up to room temperature (> 20 °C)!

### 1. Preparing the test system

- 1.1 Insert RFID card
- 1.2 Place ERS cuvette into test kit rack
- 1.3 Place ERS cap into test kit rack
- 1.4 Press measurement button and enter all required information using the touch screen

### 2. Sample preparation

- 2.1 Choose sample type: ensure that "*Pipette20µ*" is selected when performing a sample measurement. When performing a control measurement please make sure "*CtrlPlasma*" is selected.
- 2.2 Pipette 20µl sample material from centrifuged sample tube

### 3. Test processing

- 3.1 Dispense 20µl sample INTO THE LIQUID in the ERS cuvette
- 3.2 Firmly apply ERS cap on ERS cuvette
- 3.3 Place ERS cartridge into the solo analyser
- 3.4 Start automatic sample processing by pressing the  start button on the solo analyser