

The ART of Lyme Testing



Accurate, objective and automated results in **as few as 3 minutes**



Excellent performance compared to Western Blot




Flexible, dual mode testing for a variety of laboratory environments



Integrated data management automatically stores test and user history




User and patient ID captured with **onboard barcode scanner**









Sofia 2 Development Modes

WALK AWAY Mode – Walk away and multitask	READ NOW Mode – Batch multiple samples per hour

Sofia 2 Lyme+ FIA

- Differentiated IgM and IgG results with a single test
- Detection of Borrelia burgdorferi, Borrelia afzelii, and Borrelia garinii
- Objective and automated results with no equivocal zone
- Room temperature storage
- 24-month shelf life from date of manufacture

Agreement with Predicate Assay and Western Blot*

Result Type	lgM	lgG
Sensitivity	92.3%	96.4%
Specificity	90.1%	95.2%

*Refer to the Package Insert for additional performance claims.



SL2604 | SS2030900EN00 (05/19)