

# **TCM**5 FLEX MONITOR **In good hands**

Continuous and non-invasive transcutaneous monitoring of critically ill patients

## **Careful monitoring** of your most critical patients



## A non-invasive way to monitor oxygenation $(tcpO_2)$ and ventilation (tcpCO<sub>2</sub>) status in neonatal, pediatric, and adult patients.

Transcutaneous monitoring with the TCM5 FLEX monitor complements arterial and capillary blood sampling. It can help you schedule when to take blood samples and ensure careful monitoring of your most fragile patients.<sup>1,2,3,4</sup>

By being able to view continuous data instead of having to rely on periodic blood gas analysis, I feel fully confident in my ability to be precise and accurate so I can provide a higher level of care."

Dr. Marco Bardessono, Pulmonologist,

#### Focus on patient care

With its advanced features, the TCM5 FLEX helps support your workflow so you can focus on patient care.

- Easy workflow



Continuous non-invasive monitoring helps keep you in control. The TCM5 FLEX monitor gives you a real-time view of changes when they happen.

- changes in patient status

#### **Reliable sensors for accurate results**

- and greater pull strength

A.O.U. Citta della Salute e della Scienza hospital in Torino, Italy

• System defined limits of site time vs. temperature

• Built-in tutorials and clear how-to messages

Color touchscreen for a clear view of measured parameters

• Connectivity to major patient monitoring systems

#### Make informed decisions quickly and with confidence

Color touchscreen displays parameters continually

• Real-time trending measurements enable immediate clinical intervention to critical

· Clearly visible alarm bar on top of the monitor

• Surgical steel protected sensors help ensure accurate results

• A variety of fixation options support patient comfort

• Automatic calibration of the sensor reduces preparation time and improves workflow

• New for 2022: Improved sensor head design that with better protection of the cable



### Whatever comes next, we make sure life comes first

Radiometer products and solutions are used in hospitals, clinics, and laboratories in over 130 countries, providing information on critical parameters in acute care diagnostics. Through connected solutions, expert knowledge, and trusted partnership, we help health care professionals make diagnostic decisions to improve patient care.

3. Lundstrøm K. E. Continuous Blood Gas Monitoring in Neonates and Infants. neonatal INTENSIVE CARE 2003; 16, 4. 4. Janssens JP. Transcutaneous Pco2 to Monitor Noninvasive Mechanical Ventilation in Adults. Chest 1998; 113: 768-73.

Contact us today. Visit our website to arrange for a demo.