



Parameter	Comprehensive Plus	Pre-Surgical	Electrolytes	Liver	Equine / Large Animal	Equine Plus / Large Animal	Kidney	Bovine Monitoring	Kidney Monitoring	Avian Reptilian	General Health	Glucose
ALB	X	X		X	X	X	X			X	X	
ALP	X	X		X	X	X					X	
ALT	X	X		X							X	
AMY	X						X					
AST				X	X	X		X		X	X	
BUN	X	X			X	X	X		X		X	
CA	X		X		X	X	X	X		X	X	
TC				X							X	
CK					X	X		X		X		
CL			X									
TCO ₂			X				X					
CREA	X	X			X	X	X		X		X	
GGT				X	X	X					X	
GLU	X	X			X	X	X			X	X	X
K	X		X			X				X	X	
LDH					X	X						
NA	X		X			X				X	X	
MG			X			X		X				
PHOS	X		X			X	X	X	X	X	X	
TB	X			X	X	X					X	
TBA				X						X		
TG					X	X						
TP	X	X		X	X	X				X	X	
URIC ACID										X		
BUN/CREA*	X	X			X	X	X		X		X	
GLOB*	X	X		X	X	X				X	X	
A/G*	X	X		X	X	X				X	X	

*berechnete Parameter